Digital Microscope VHX-6000 Series

KEYENCE



KEYENCE

VHX



Built-in brightfield, darkfield, and mixed lighting to view any sample

Combines the lighting techniques of a metallurgical microscope (brightfield) and stereoscope (darkfield) so that optimal lighting can be achieved for any sample. Exact lighting conditions can be saved and recalled later for fast, reproducible imaging.



BRIGHTFIELD

The surface of the package area reflects too much light, and it cannot be seen clearly.



DARKFIELD

The light is scattered and not enough light is received by the lens to observe the sensor.



MIXED LIGHTING

Combining brightfield and darkfield lighting allows for clear observation of the entire sample.

Easy observation, even at an angle







The compact design of the dual objective zoom lens along with the flexibility of the free-angle stand means that you can easily view your sample from any angle, at any magnification.

Model		VH-ZST							
Magnification ¹		20×	100×	200×	500×	1000×	2000×		
	H (Horizontal)	15.24 0.60"	3.05 0.12"	1.52 0.0598"	0.61 0.0240"	0.30 0.0118"	0.15 0.0059"		
Field-of-view (mm inch)	V (Vertical)	11.4 0.45"	2.28 0.0898"	1.14 0.0449"	0.46 0.0181"	0.23 0.0091"	0.11 0.0043"		
	D (Diagonal)	19.05 0.75"	3.81 0.15"	1.91 0.0752"	0.76 0.0299"	0.38 0.0150"	0.19 0.0075"		
Working distance (mm	inch)	15,0,59							

*1 Magnification using a 1/2-inch CCD camera and 15-inch monitor.

Industry-Leading High-Resolution RZ Lenses



TRIPLE'R-compliant lenses are fitted with Automatic Lens/Zoom Recognition units.

High-Performance Low-Range Zoom Lens **VH-ZOOR/ZOOT**



20 200

1001000

250 2500

500 5000

Macro zoom lens

With a magnification range from 0.1× to 50×, a target can be viewed from its entirety down to more in-depth observation. This macro lens excels in workability and high performance with click-style magnification adjustment, an aperture mechanism, and a viewing distance of 95 mm 3.74° or more.

Model				VH	-Z00R/Z0	IOT		
Magnif	ication*	0.1×	0.5×	1×	5×	10×	30×	50×
Ma (Horizontal	3200 125.98"	640 25.20"	320 12.60*	61 2.40"	30.5 1.20"	10.2 0.40"	6.1 0.24
Field of view (mm inch)	Vertical	2400 94.49"	480 18.90"	240 9.45"	45.5 1.79"	22.8 0.90"	7.6 0.30"	4.6 0.18
Eiel (m	Diagonal	4000 157.48"	800 31.50"	400 15.75*	76.2 3.00"	38.1 1.50"	12.7 0.50"	7.6 0.30'
Working distance (mm inch)		Approx. 7700 303.15"	Approx. 1500 59.06"	Approx. 720 28.35*		9 3.1	5 74"	

* Magnification on a 15-inch monitor with a 1/2-inch CCD camera.

Ultra-Small, High-Performance Zoom Lens VH-Z20R/Z20T

High-resolution imaging with a large depth-of-field

The VH-Z20R/Z20T offers high-resolution observation at general purpose magnifications of 20× - 200x. This lens has been designed to optimize both depth-of-field and resolution and can also be used in handheld mode.

High-performance lens with long working distance This innovative lens was developed to satisfy the need for high-resolution, long working distance, and large depth-of-field.

Model	VH-Z20R/Z20T						
Magnif	ication*1	20×	30×	50×	100×	150×	200×
) ew	Horizontal	15.24 0.60"	10.16 0.40"	6.10 0.24"	3.05 0.12	2.03 0.0799"	1.52 0.0598"
Field of view (mm inch)	Vertical	11.40 0.45"	7.60 0.30"	4.56 0,18"	2.28 0.0898"	1.52 0.0598"	1.14 0.0449"
	Diagonal	19.05 0.75"	12.70 0.50"	7.62 0.30"	3.81 0.15"	2.54 0.10"	1.91 0.0752"
Depth of field (mm inch)*2		34 1,34"	15.5 0.61"	6.0 0.24	1.6 0.0630"	0.74 0.0291"	0.44 0.0173"
Workin (mm in	g distance ch)				5.5 00"		

*1 Magnification on a 15-inch monitor with a 1/2-inch CCD camera.
*2 This value is for when the priority is set to depth of field. The depth of field changes depending on the set aperture.

Wide-Range Zoom Lens VH-Z100R/Z100T

Darkfield

Model				VH-Z100	R/Z100T		
Magnif	ication*1	100×	200×	300×	500×	700×	1000×
ew (Horizontal	3.05 0.12"	1.53 0.0602"	1.02 0.0402"	0.61 0.0240"	0.44 0.0173"	0.30 0.0118"
Field of view (mm inch)	Vertical	2.28 0.0898"	1.14 0.0449"	0.76 0.0299"	0.46 0.0181"	0.33 0.0130"	0.23 0.0091*
	Diagonal	3.81 0.15"	1.90 0.0748"	1.27 0.0500"	0.76 0.0299"	0.54 0.0213"	0.38 0.0150*
Working distance (mm inch) 25 0.98" (20 0.79"*2)							

*1 Magnification on a 15-inch monitor with a 1/2-inch CCD camera. *2 With triple illumination adapter attached

Brightfield Polarization

Dual-Light High-Magnification Zoom Lens VH-Z250R/Z250T

Observe with both brightfield and darkfield at high-magnification

Easily switch between ring and coaxial illumination with just the touch of a button. View objects at up to $2500\times$ magnification while still maintaining a 6.5 mm 0.26° working distance.

Brightfield Darkfield

Model			.22 1.02 0.61 0.31 0.2 0.15 0.12 (480" 0.0402" 0.0240" 0.0122" 0.0079" 0.0059" 0.0047 92 0.76 0.46 0.23 0.15 0.11 0.09							
Magnif	ication*	250×	300×	500×	1000×	1500×	2000×	2500×		
iew I	Horizontal	1.22 0.0480"						0.12 0.0047"		
Field of view (mm inch)	Vertical	0.92 0.0362"						0.09 0.0035"		
	Diagonal	1.52 0.0598"	1.27 0.0500"	0.76 0.0299"	0.38 0.0150'	0.25 0.0098"	0.19 0.0075"	0.15 0.0059"		
Working distance (mm inch)			6.5 0.26"							
* * * * * * * *	the set of									

Magnification on a 15-inch monitor with a 1/2-inch CCD camera.

High-Resolution Zoom Lens VH-Z500R/Z500T

Our highest magnification/ resolution zoom lens

This zoom lens incorporates high-quality fluorite optics to provide the highest resolution in its class. With an N.A. of 0.82, achieve up to $5000 \times$ magnification with a 4.4 mm 0.17° working distance.

Model VH-Z500R/Z500T						
Magnification*		500×	1000×	2000×	3000×	5000×
ъĉ	Horizontal	610	305	152	102	61
Field of view (µm)	Vertical	457	229	114	76	46
Kie Fi	Diagonal	762	381	191	127	76
Working distance (mm inch)				4.4 0.17 "		

* Magnification on a 15-inch monitor with a 1/2-inch CCD camera.









SPECIFICATIONS

Basic functions: Controller

Model			VHX-6000	VHX-950F
	Image senso	r	1/1.8-inch CMOS image sensor Virtual pixels: 1600 (H) × 1200 (V)	1/1.8-inch CMOS image sensor Virtual pixels: 1600 (H) × 1200 (V)
	Scanning sys	stem	Progressive	Progressive
	Frame rate		50 F/s (max.)	50 F/s (max.)
	N	orma	1600 (H) × 1200 (V), Approx. 1000 TV lines	1600 (H) × 1200 (V), Approx. 1000 TV lines
	-	CMOS ^{*1*3}	1600 (H) × 1200 (V), Approx. 1200 TV lines or more (2 million pixels × 3CMOS mode, Excellent color reproducibility)	
	Image resolution	ligh resolution*3	3200 (H) × 2400 (V), Approx. 1600 TV lines	Natawiahla
2010010	S	uper Fine ^{*3}	4800 (H) × 3600 (V), Approx. 2000 TV lines or more	Not available
Camera	s	uper Fine × 3CMOS ^{-2'3}	4800 (H) × 3600 (V), Approx. 2000 TV lines or more (18 million pixels × 3CMOS mode, Excellent Color reproducibility)	
	High Dynami	c Range	16-bit intensity range through RGB data from each pixel	Not available
	Gain		Auto, Manual, Preset	Auto, Manual, Preset
	Electronic sh	utter	Auto, Manual, 1/60, 1/120, 1/250, 1/500, 1/1000, 1/2000, 1/5000, 1/9000, 1/19000	Auto, Manual, 1/60, 1/120, 1/250, 1/500, 1/1000, 1/2000, 1/5000, 1/9000, 1/1900
	Supercharge	shutter	0.02 to 4 s	0.02 to 4 s
	White balance	e	Push set, Auto, Manual, Preset (2700K, 3200K, 5600K, 9000K)	Push set, Auto, Manual, Preset (2700K, 3200K, 5600K, 9000K)
	Back-focus a	djustment	Not required	Not required
	Size	·	Color LCD (IPS) 23"	Color LCD (IPS) 23"
	Screen size		509 184 (H) × 286 416 (V) mm 20.05"(H) × 11.28"(V)	509.184 (H) × 286.416 (V) mm 20.05"(H) × 11.28"(V)
	Pixel pitch		0.2652 (H) × 0.2652 (V) mm 0.01"(H) × 0.01"(V)	0.2652 (H) × 0.2652 (V) mm 0.01"(H) × 0.01"(V)
	Number of p	ixels	1920 (H) × 1080 (V) (FHD)	1920 (H) × 1080 (V) (FHD)
LCD monitor*5	Display colors		Approx. 16770000 colors*4	Approx. 16770000 colors*4
	Brightness		300 cd/m ² (Center 1 Point, typical)	300 cd/m ² (Center 1 Point, typical)
	Contrast ratio		1000:1 (typical)	1000:1 (typical)
	Field of view		±89° (typ., horizontal), ±89° (typ., vertical)	±89° (typ., horizontal), ±89° (typ., vertical)
	Unit		DVD-ROM super-multi drive unit	DVD-ROM super-multi drive unit
CD-R/CD-RW/DVD	Applicable disk		CD-R/CD-RW/DVD±R/DVD±R DL/DVD±RW/DVD-RAM	CD-R/CD-RW/DVD±R/DVD±R DL/DVD±RW/DVD-RAM
drive unit	Storage capa	citv	8.7 GB (when using DVD±R DL)	8.7 GB (when using DVD±R DL)
Hard disk drive unit	Storage capa	city	500 GB (including 165 GB reserved area) Approx. 1680000 images (when a 2-million-pixel image is compressed) to approx. 55000 images (when a 2-million-pixel image is not compressed)	500 GB (including 165 GB reserved area) Approx. 1680000 images (when a 2-million-pixel image is compressed) to approx, 55000 images (when a 2-million-pixel image is not compressed)
mage format			JPEG (with compression), TIFF (without compression)	JPEG (with compression), TIFF (without compression)
Observable image siz	e		20000 (H) × 20000 (V) pixels (when stitched)	1600 (H) × 1200 (V) pixels
obool vabio iniago biz	Lamp		High-brightness LED	High-brightness LED
Light source	Lamp life		40000 hours (reference)	40000 hours (reference)
-igitt oouroo	Color temper	ature	5700K (typical)	5700K (typical)
	Output meth		DVI-I (1920 × 1080 pixels)	DVI-I (1920 × 1080 pixels)
Video output	Scanning	Special LCD monitor	66 kHz (H), 60 Hz (V)	66 kHz (H), 60 Hz (V)
naco output	frequency	External monitor	66 kHz (H), 60 Hz (V)	66 kHz (H), 60 Hz (V)
	Mouse input		USB mouse supported	USB mouse supported
Input	Keyboard ing		USB keyboard supported	USB keyboard supported
mput	External rem		Pause/Recording, Non-voltage input (Contact/Noncontact)	Pause/Recording, Non-voltage input (Contact/Noncontact)
	LAN		RJ-45 (10BASE-T/100BASE-TX/1000BASE-T)	RJ-45 (10BASE-T/100BASE-TX/1000BASE-T)
nterfaces	USB 2.0 Seri	es A	6 types	6 types
mornaooo	USB 3.0 Seri		2 types	2 types
	Power voltag		100 to 240 VAC, 50/60 Hz	100 to 240 VAC, 50/60 Hz
Power supply	Power consu		280 VA	280 VA
Environmenta	Ambient tem		+5 to 40°C 41 to 104°F	+5 to 40°C 41 to 104°F
resistance	Relative hum		35 to 80% RH (No condensation)	35 to 80% RH (No condensation)
osistanos	Controller	nuny		
Weight	Camera unit		Approx. 12.5 kg Approx. 1.2 kg (VHX-6100/6020)	Approx. 12.5 kg Approx. 1.2 kg (VHX-6020)
weight	Camera unit		Approx. 1.2 kg (VHX-6100/6020) Approx. 0.5 kg	Approx. 1.2 kg (VHX-6020) Approx. 0.5 kg
	UCONSOLE		Approx. 0.5 kg $550 (W) \sim 470 (H) \sim 200 (D) 21.65"(W) \sim 18.50"(H) \sim 7.87"(D) (when stored)$	Approx. 0.5 kg 550 (W) × 470 (H) × 200 (D) 21 65°(W) × 18 50°(H) × 7.87°(D) (when stored)

Dimensions (Excluding the projected areas) 550 (W) × 470 (H) × 200 (D) 21.65° (W) × 18.50° (H) × 7.87° (D) (when stored) 550 (W) × 470 (H) × 200 (D) 21.65° (W) × 18.50° (H) × 7.87° (D) (when stored) Т *1. Compared with the standard mode, resolution and color reproducibility are improved.
 *2. Compared with the high-resolution HD mode, color reproducibility is improved.
 *3. Supported only with the VHX-6100 multi-scan camera.
 *4. Approx. 16.77 million colors can be reproduced with the FRC processing of the display controller.
 *5. The LCD monitor provided with the VHX Series is based on advanced technology. Rarely, an unlit pixel (black spot) or lit pixel (bright spot) may exist on the monitor screen. However, this is not an indication of the LCD monitor being defective.

Basic functions: Stage

		VHX-S660E	VHX-S650E	VHX-S600E	VH-S300	VH-S30F/S30B
	XY stage: Motorized/Manual	Motorized	Motorized	Manual	Manual	Manual
	XY-motorized stage motor	2-phase stepping motor	2-phase stepping motor	-	-	-
	XY-motorized stage resolution	1 µm (typ.)	1 µm (typ.)	-	-	-
	XY-motorized stage movement speed	20 mm 0.79"/sec (max.)	10 mm 0.39"/sec (max.)	-	-	-
$XY\theta$ stage	XY-stage movement range	±50 mm ±1.97"	±20 mm ±0.79"	±35 mm ±1.38"	±35 mm ±1.38"	X: ±37.5 mm ±1.48", Y: ±25 mm ±0.98"
Aro olago	θ rotation angle	-	±90°	360°	360°	360°
	XY $ heta$ stage size	Top surface: 233 × 185 mm 9.17" × 7.28" (Center disc: ø168 ø6.61")	Top surface: 171 × 168 mm 6.73" × 6.61" (Center disc: ø100 ø3.94")	Top surface: 198 × 150 mm 7.80" × 5.91" (Center disc: ø136 ø5.35")	Top surface: 190 × 150 mm 7.48" × 5.91"	Top surface: 180 × 136 mm 7.09" × 5.35"
	Transmitted lighting	Provided	Provided	Provided	-	-
	Z stage: Motorized/Manual	Motorized	Motorized	Motorized	Manual	Manual
	Z-motorized stage motor	5-phase stepping motor	5-phase stepping motor	5-phase stepping motor	-	-
Z stage	Z-motorized stage resolution	0.1 µm (typ.)	0.1 µm (typ.)	0.1 µm (typ.)	-	-
	Z-motorized stage travel speed	17 mm 0.67"/sec (max)	17 mm 0.67"/sec (max)	17 mm 0.67"/sec (max)	-	-
	Z-stage movement range	49 mm 1.93"	49 mm 1.93"	49 mm 1.93"	53 mm 2.09"	28 mm 1.10"
Ratings	Power voltage	100 to 240 VAC, 50/60 Hz	100 to 240 VAC, 50/60 Hz	-	-	-
naunys	Power consumption	60 VA	60 VA	50 VA	-	-
Environmental	Ambient temperature	+5 to 40°C 41 to 104°F	+5 to 40°C 41 to 104°F	+5 to 40°C 41 to 104°F	-	-
resistance	Relative humidity	35 to 80% RH (No condensation)	35 to 80% RH (No condensation)	35 to 80% RH (No condensation)	-	-
Weight		Approx. 20 kg	Approx. 18 kg	Approx 17.2 kg	Approx. 17.4 kg	Approx. 12 kg
Load capacity		5 kg	5 kg	5 kg	1 kg	1 kg

Other functions

Model		VHX-6000		VHX-950F	Console compatibil
	Auto-focus function	Available	✓ ✓	Available	1 V
	Focus amount display function	Available		Available	
Observation functions	Name of the second se	1			
	Multi-lighting function	Available	compatibility Unresour Available Available Available Available Sembination lighting) - Not available Available - Available Available - Available - Available - - Not available - - Available - <td></td>		
	Camera-shake correcting function	Available	1	Available	1
	Full-screen display	Available	1	Available	1
Diaplay function	Screen segmentation function	Vertical, horizontal, 2 × 2 split, 3 × 3 split, linked display function		Vertical, horizontal, 2 × 2 split, 3 × 3 split, linked display function	1
Jispilay function	Real-time digital zoom	1.0× to 10.0×		1.0× to 10.0×	1
	Comment display function	Available		Available	
	Glare removal function	Available	1	Available	1
	Ring removal function		1	Not available	
mage qua l ity			1		
mprovement	High-resolution HDR function		1		
Focus answer (bing surface frame) Available Image: Available Image: Available Image: Available Commentation function Gall, partal, later, darield, later, daried, later, darield, later, daried, later, darield, later, daried					
	Fine-shot function	Available	1	Available	1
	2D image stitching	Available	1	Not available	
Stitching function	3D image stitching	Available	1	Not available	
	Navigation function	Available		Not avai lable	
	Real-time depth composition	Available	1	Not avai lable	
	Quick composition & 3D function	Available	1	Available	1
3D function					
	3D display function	Available		Available	
	3D shape correction function	Available (Slope/Sphere/Cylinder)			
	3D comparison function	Available (Combination/Comparison/Difference display mode)		Available (Combination/Comparison/Difference display mode)	
	Report output (Excel/Word)				
	Capture setting reproduction function	Available		Available	
Recording function	Timer capture function	Available Available nction Available Available Available Available Available S0 frames/sec. max. 30 frames/sec. max. (Image size: 1600 × 1200, 800 × 600, 640 × 480) (Image size: 1600 × 1200, 800 × 600, 640 × 480) c. Variously available Variously available Available Variously available Variously available			
	Video recording/playback function				
	Distance, angle, radius, area, etc.	Variously available		Variously available	1
	Automatic edge detection	Available		Available	
	Scale display	Variously available	1	Variously available	1
7					
	Automatic area measurement	Available		Not available	
	Grain size analysis	Available		Not available	
Measuring functions	Contamination analysis	Available		Not available	
	One-click measurement	Available		Not available	
	Auto measurement	Available		Not available	
	function (TRIPLE'R)				
	<u> </u>				
3D measurement					
unction (VHX-H4M/	×				
					ļ
,	-				
measurement system (VHX-H2M2 optional					<u> </u>
function)					
			*		- *
Itilities					
Junios					
	Network connection function	Available (Communication software / File sharing / FTP)		Available (Communication software / File sharing / FTP)	
			+		
			-		
	So image playback software for the PC	, , , , , , , , , , , , , , , , , , , ,		heproduces 3D images saved on the VHX using the PC.	
PC software (Available free of		already been captured.		Not available	
/HX-S15 optional unction) Wanual XY neasurement system (VHX-H2M2 optional unction) Jtilities 2C software	HDR playback/measurement/stitched image playback software			Performs measurement using the PC.	
	One-click measurement compilation		1	N 1 1 1 1	1

VHX-6000/950F (Module details)

	Video recording module	Allows recording/playback of videos.			
	High-quality depth composition module	Captures multiple images focused on different heights from which a single image is composed			
	Area measurement module	easures an area of a 2D image.			
Module	Timer capture module	Captures images automatically at specified time intervals.			
	Screen splitting module	Displays vertical, horizontal, or four-division split screens.			
	Comment input module	Allows inputting and displaying of comments such as characters and markers on an image.			
	Image touch-up module	Provides image processing functions for modifying images to make observation easier.			



www.keyence.com



SAFETY INFORMATION

Please read the instruction manual carefully in order to safely operate any KEYENCE product.

CONTACT YOUR NEAREST OFFICE FOR RELEASE STATUS

KEYENCE CORPO Head Office 500			a, IL 60143, U.S.A.	PHONE: +1-201	-930-0100 FAX: +1 [,]	-855-539-0123	E-mail: keyence@ke	yence.com	
AL Birmingham AR Little Rock AZ Phoenix CA San Francisco	CA San Jose CA Cupertino CA Los Angeles CA Irvine	CO Denver FL Tampa GA Atlanta IA Iowa	IL Chicago IN Indianapolis KY Louisville MA Boston	MI Detroit MI Grand Rapids MN Minneapolis MO Kansas City	MO St. Louis NJ Elmwood Park NY Rochester NC Charlotte	NC Raleigh OH Cincinnati OH Cleveland OR Portland	PA PhiladelphiaPA PittsburghSC GreenvilleTN Knoxville	TN Nashville TX Austin TX Dallas WA Seattle	WI Milwaukee
	IE: +1-905-366-765		66-1122 E-mail: keye 94-3206 Windsor PH		e.com 55 FAX: +1-905-366-	.1122	KEYENCE MEXICO PHONE: +52-55-8850 E-mail: keyencemexio	0-0100 FAX: +52	-81-8220-9097

The information in this publication is based on KEYENCE's internal research/evaluation at the time of release and is subject to change without notice. Company and product names mentioned in this catalog are either trademarks or registered trademarks of their respective companies. The specifications are expressed in metric units. The English units have been converted from the original metric units. Copyright © 2017 KEYENCE CORPORATION. All rights reserved.

KA1-1017