



**COVALENT
METROLOGY**

CM000030169

Covalent Metrology, Max Junda

CdS Bandgap Mapping by Spectroscopic Ellipsometry

Authors Max Junda, PhD

Reviewer Avery Green, PhD

Published Date September 1, 2023

Goal: *Measurement of the spatial uniformity of the bandgap of a CdS thin film both before and after annealing.*

Methods: Mapping spectroscopic ellipsometry is collected in a 36-point grid across the full sample area before and after annealing.


Results Summary:

- The average bandgap energy is found to increase after annealing with the following averages:
 - Before Annealing: **2.23 eV**
 - After Annealing: **2.37 eV**
- The bandgap spatial uniformity is decreased following annealing.
- The increased surface roughness following annealing suggests grain growth as a result of the annealing process.
- Raw results data is provided alongside this report in spreadsheet format.

Measurement Results


Before Anneal

Bandgap (eV) vs. Position




After Anneal

Bandgap (eV) vs. Position




- The mean squared error function (MSE) is a relative measure of how well the model fits the measured data.
- It quantifies any mismatch between the model-generated and measured spectra such that lower MSE values correspond to better model fits.



- The surface roughness of the film approximately doubles as a result of annealing.
- Generally, annealing is known to improve the crystallographic order in a material and promote the growth of larger grains.
- Increased surface roughness is consistent with larger grain size since larger individual grains protruding from the surface more.
- Additional characterization (such as XRD) would be required to confirm this grain size hypothesis with certainty.


Before Anneal


Roughness (nm) vs. Position



After Anneal

Roughness (nm) vs. Position





Before Anneal


Statistic	Value
Average	2.2304
Std. Dev.	0.0026
Minimum	2.2244
Maximum	2.2339
Range	0.0095
# of Points	36

After Anneal

Statistic	Value
Average	2.3724
Std. Dev.	0.0042
Minimum	2.3653
Maximum	2.3786
Range	0.0134
# of Points	36


Before Anneal

Spectroscopic Data at X=2.25, Y=0



After Anneal

Spectroscopic Data at X=-0.75, Y=1.5



+ Layer # 2 = Roughness Roughness Thickness = 9.31 nm (fit)
+ Layer # 1 = CdS CdS Thickness = 96.55 nm (fit)
Substrate = Glass

- For both the pre- and post-anneal measurements, the three-layer model shown above is fit to the measured ellipsometric spectra at each mapping location.
- The optical properties of the CdS film are freely fit at each location and are described as the sum of three Tauc-Lorentz and one Sellmeier oscillator.
- Tauc-Lorentz oscillators are typically used to describe amorphous or disordered semiconductors and feature a variable bandgap parameter which controls the photon energy where the optical absorption goes to zero – the value of this parameter is what is mapped to determine spatial bandgap variation.
- The surface roughness is described as a Bruggeman effective medium approximation consisting of a 50/50% mix of CdS material and void.

Appendix 1

J. A. Woollam Co., Inc. Model RC2-DI Spectroscopic Ellipsometer

- Dual rotating compensator configuration multichannel spectroscopic ellipsometer
- Full spectrum data acquisition in:
 - Minimum: 0.3 s
 - Typical: 10 – 20s
- Automated mapping up to 300 mm substrates
- 190 to 1690 nm spectral range
- CompleteEase modeling software
- Variable angle transmission stage
- Measurement beam diameter:
 - Normal: 5 mm
 - Focused: 300 μm



Order


[New Quote](#)



Order Status


<u>Order Number</u> ▾	<u>Accepted Date</u> ⁱ ▾	<u>Samples Received</u> ▾	<u>Status</u> ⁱ ▾	
CM000023633 1 samples, 2 techniques	2022-06-24	1	● Waiting for samples ⁱ	Repeat
CM000023632 1 samples, 2 techniques	2022-06-24	1	● Sample Received ⁱ	Repeat
CM000023631 1 samples, 2 techniques	2022-06-24	1	● Sample Received ⁱ	Repeat
CM000021056 1 samples, 2 techniques	2022-01-25	1	● Executing ⁱ	Repeat
CM000020337 0 samples, 4 techniques	NA	0	● Waiting for samples ⁱ	Repeat
CM000020336 1 samples, 4 techniques	2022-01-03	1	● Waiting for samples ⁱ	Repeat

The Solution Portal can be accessed from your Order Portal dashboard

Watch the Demo

- **Interact with Images, PPTX Files, and EDX/EDS Data**
- **Download & Save** analysis, measurements and results
- **Pre-built Analytical Tools** accelerate analytical routines and reporting
 - Line Edge Roughness (LER): interfacial roughness
 - Area Measurement: particle and grain size analysis
 - Line Profile: thin film thickness and multilayer measurements

Start Exploring at:
solution.covalentmetrology.com



The screenshot displays the software interface with several callout boxes highlighting key features:

- Share important notes to your team in the Cloud!**: A callout box pointing to the 'Notes' field in the top toolbar.
- Copy your analysis into reports or presentations**: A callout box pointing to the 'Copy' icon in the top toolbar.
- Make your own measurements and annotate your images**: A callout box pointing to the measurement tool being used on the image.
- Access & Analyze your Images as they are collected**: A callout box pointing to the file gallery on the left.
- Favorite Results to find your most important data fast**: A callout box pointing to the star icon in the file gallery.
- Dynamic Scalebar updates as you crop and zoom your image**: A callout box pointing to the '20 nm' scalebar at the bottom right.

The interface includes a top menu bar with 'File', 'Edit', 'Annotation', and 'Tools'. Below it are toolbars for text, measurement, and notes. The main area shows a file gallery on the left and a large image on the right with measurement lines and a scalebar.



**COVALENT
METROLOGY**

Thank You.

covalentmetrology.com